



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fuel Injector Cleaner And Tester Four Injectors

Product Code :
LBNY-0005-12100022



Description :

The device simulates various working conditions of the engine, cleaning and detecting the fuel injection nozzle of various automobile, and can also carry out the cleaning of the automobile fuel injection nozzle and the fuel supply system. The equipment adopts ultrasonic cleaning technology and microcomputer oil pressure closed loop control cleaning detection technology.

Features:

- Ultrasonic Cleaning
- On-vehicle Cleaning
- Auto cleaning test function
- Uniformity & Spray ability Test
- Leakage Test
- Direct Flush
- Reverse Flush

Technical Specification :

Simulated RPM range: 0-9990rpm, step: 10rpm
Power Supply: AC 220v~ 50HZ
Pressure: 0-6bar
PM range: 0.1- 20ms, Step: 0.1ms
Input Power: 200W (minimum)



Laboratory Instrument India