



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Inverter Air Conditioner Training Equipment

**Product Code :**  
LBNY-0005-12100024



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Inverter Air Conditioner Training Equipment can be used to train students to master the theoretical knowledge and refrigeration professional practice skills of inverter air conditioning system.

**Technical Specification :**

Dimension: about 950\*800\*1750mm  
Power supply: AC220V, 50HZ, °1.5KW.



**Laboratory Instrument India**