



Email : sales@laboratoryinstrumentindia.com

Product Name :
PLC Programmable Controller Training Platform

Product Code :
LBNY-0005-12100027



Description :

Modular structure, hanging box is easy for doing different experiment, easy to organize various training project, expand and upgrade are convenient.

Main Functions:

It has a safety protection function.

It can finish PLC commonly used programming instructions learning and training, PLC control car linear motion project training, PLC typical application project training, PLC control three-phase asynchronous motor running project training.

Technical Specification :

Power supply: three-phase AC380V, power Size: 1700 mm x 700 mm x 1500 mm



Laboratory Instrument India