



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Inverter Air Conditioner Circuit Maintenance Training Device

**Product Code :**  
LBNY-0005-12100028



**Description :**

Inverter Air Conditioner Circuit Maintenance Training Device can be used to train students to master the theoretical knowledge and practice skills of inverter air conditioning system.

**Technical Specification :**

Input power: AC220V  $\pm$  10%, 50Hz.

Total faults: 18 .

Power supply: 220V  $\pm$  10%, 50Hz.

Output: 1 channel AC 220V output.

a. 1 channel regulated output DC 5V

b. 1-channel adjustable output 0 ~ 5V DC

Dimensions: approximately 1400 x 700 x 1700mm ( L x W x H ) .

Training Bench adopt industrial profiles, high-quality, light weight, with four universal wheel brakes, it can be easily moved.



**Laboratory Instrument India**

