



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air Conditioner Installation And Maintenance Training Device

Product Code :
LBNY-0005-12100033



Description :

Air Conditioner Installation And Maintenance Training Device can be used to train students to master the theoretical knowledge and practice skills of design, installation, commissioning and maintenance of air conditioning system.

Technical Specification :

Structure: The entire system is mounted on the mobile aluminum profile frame with universal wheel.
Rated power: 1.5KW
Input power: AC220V, 50HZ
Dimension: about 950*800*1750mm



Laboratory Instrument India