



Email : sales@laboratoryinstrumentindia.com

Product Name :
Small Freezer Training Equipment

Product Code :
LBNY-0005-12100035



Description :

The students can learn the cold storage unit design, load calculation, electrical design, warehouse design, cold storage management, cost estimates and other aspects skills. It can simulate the control principle and workflow of real small freezer. Small Freezer Training Equipment includes three main parts of electrical control systems, refrigeration units, and real storage room.

Technical Specification :

Safety protection: overvoltage, over current, overload, leakage and grounding.
Compressor power: 2.5HP
Refrigerant: R22
Power supply: AC380V \pm 10%, 50Hz
Dimension: approximately 6000mm*2400mm*2500mm



Laboratory Instrument India