



Email : sales@laboratoryinstrumentindia.com

Product Name :
Inverter Air Conditioner Training Equipment

Product Code :
LBNY-0005-12100037



www.laboratoryinstrumentindia.com

Description :

Inverter Air Conditioner Training Equipment can be used to train students to master the theoretical knowledge and refrigeration professional practice skills of inverter air conditioning system.

Technical Specification :

Dimension: about 950*800*1750mm
Power supply: AC220V, 50HZ, 1.5KW



Laboratory Instrument India