



Email : sales@laboratoryinstrumentindia.com

Product Name :
Analog Circuit Training Kit

Product Code :
LBNY-0005-12100038



Description :

All the experiment module are located on the circuit board, the trainees can build the experiment circuit by connecting the different module via locking jack leads.

Experiments:

- Differential amplifier circuit
- Proportional summation circuit
- Integral and differential circuit
- Negative feedback amplifier circuit
- Emitter follower
- Waveform generating circuit
- Single-stage amplifier circuit
- Two-stage amplifier circuit
- Active filter.

Technical Specification :

Analog Circuit Training Kit

Laboratory Instrument India

