



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Refrigerator Repair And Training Equipment Direct Cooling

**Product Code :**  
LBNY-0005-12100039



**Description :**

Refrigerator Repair And Training Equipment Direct Cooling can be used to train students to master the theoretical knowledge and practice skills of direct cooled refrigerator system.

The control panel wiring connection diagram panel with test points in the circuit diagram, it is easy for students to measure and debugging.

Structure And Compositions:

Mobile aluminum profile frame with control panel:

Control panel wiring connection diagram panel,

Refrigerators refrigeration schematics panel,

Fault setting panel.

**Technical Specification :**

Rated power: 1.0KW

Power supply: 220V AC, 50Hz;

Dimensions: about 950mm\*800mm\*1600mm



**Laboratory Instrument India**

