



Email : sales@laboratoryinstrumentindia.com

Product Name :
Analog Multimeter Electronic Engineering Lab

Product Code :
LBNY-0002-620045



Description :

Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab

Technical Specification :

A typical multimeter would include basic features such as the ability to measure voltage, current, and resistance. Analog multimeters use a microammeter whose pointer moves over a scale calibrated for all the different measurements that can be made.

Features:

- Analogue Multimeter with Mirror Scale
- Inbuilt Buzzer and LED
- Transistor Test, hFE with Diode

Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab, Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab Equipments, Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab Tools, Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab Tool Kits, Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab Manufacturers, Analog Multimeter (DC 1000V, AC 750V, 1K?) Electronic Engineering Lab Suppliers from India, China, Kenya

Laboratory Instrument India

