



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Hydraulic Transmission Training Device

Product Code :
LBNY-0005-12100043



Description :

It is suitable for vocational colleges related courses such as hydraulic transmission, hydraulic and pneumatic transmission of practice teaching. Portable Hydraulic Transmission Training Device transparent hydraulic transmission experiment box is mainly for demo of the hydraulic transmission and hydraulic transmission basic loop (pressure, flow, direction control loop).

Features:

Hydraulic components, piping in a transparent hose, compact structure, reliable connection and no leakage. Hydraulic components in a transparent organic glass shell as the basic, working medium for hydraulic oil, white system at work can be clearly observed that the working process of the hydraulic components of the valve core, intuitive is strong.

Box set element component, hydraulic pump station, control motor and electric control part of box, the teacher may convenient will be live demonstration equipment into classroom teaching and the physical interpretation.

Technical Specification :

Oil pump: rated pressure 1.6 MPa, rated flow: 4 l/min
Power supply: Dc motor speed control: AC220V, 180 w, 1400 r/min
Outside size: 600 mm * 400 mm * 600 mm.

Laboratory Instrument India

