



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Basic Hydraulic Training Device

**Product Code :**  
LBNY-0005-12100047



#### **Description :**

Basic Hydraulic Training Device adopts transparent hydraulic components, profiles aluminum alloy structure, active oil joint and general electric circuit line. With the use of attached tool, it is convenient to control and do most experiments in hydraulic transmission and it can let the students understand the principles, materials and the experimental instruction, structure and working process of oil circuit and hydraulic components internal structure. It is an essential instrument for hydraulic transmission teaching.

#### **Technical Specification :**

Electrical control components: Contains:  
Power supply unit, Relay control unit, circuit protection unit, the experimental station protection unit, motion valve (normally open, normally closed) and so on.  
Input voltage: AC220V, 50HZ  
Output voltage: DC24V, AC24V  
Dimensions: 1500mm(L)\*700mm(W)\*1700mm(H).



**Laboratory Instrument India**

