



Email : sales@laboratoryinstrumentindia.com

Product Name :
Refrigerant Charging Station With Scale

Product Code :
LBNY-0005-12100048



Description :

Features:

Portable and easy to carry
Built in safety valve manifold compound and pressure gauges
Vacuum gauge hose ccl-36 (90cm) with safety goggles and oil for vacuum pump
Compressor power: 280W
Rotary vane type with two stage vacuum pump
With a digital scale (up to 100 Kg with a 0.5% accuracy)
Vacuum pump pumping speed: 60L/min
Filling speed: 800g/min
Recovery rate: 150g/min
Working tank volume: 12L
External hose: 2.0M
Refrigerant electronic weighing range: 60Kg
Refrigerant electronic weighing accuracy: + 10g
Temperature: 0 -40 degree, relative humidity: 85%
External magnetic field strength: 400mA/m
Power supply: AC220 + 10% frequency: 50/60HZ
Package size: 620X680X1160mm
Machine weight: 90 kg.

Technical Specification :

Refrigerant Charging Station With Scale



Laboratory Instrument India