



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab

**Product Code :**  
LBNY-0002-620046



**Description :**

Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab

**Technical Specification :**

Digital Multimeter (AC 500V, DC 500V, 200M?)

It's a good electronic device. Stylish appearance, easy to use.

LCD display design makes you can fast read data. It can be used to measure voltage current resistance, etc. It has stable and high reliability with an anti-drop handle.

This instrument has a backlight, user-friendly read the measure value display in a dark place.

It can be used to measure DC and AC voltage, DC current resistors, diodes, transistors and circuit continuity test.

Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab, Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab Equipments, Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab Tools, Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab Tool Kits, Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab Manufacturers, Digital Multimeter AC 500V DC 500V 200M? Electronic Engineering Lab Suppliers from India, China, Kenya

---

# Laboratory Instrument India

