



Email : sales@laboratoryinstrumentindia.com

Product Name :
Four Way Manifolds Digital

Product Code :
LBNY-0005-12100052



Description :

Pressure and vacuum mode shift function build-in 82 refrigerants according to standard refrigerants pressure and temperature relationship database.

Real-time vacuum pressure test and percent indicate

Test range Measurement of vacuum: -101Kpa~0Mpa

Positive pressure: 0Kpa~6Mpa

Low power indicate. Batteries life:30H. 9V.

Pressure units: KPa, Bar, Psi, Torr, Kg F/cm²,InHg, MmHg

Refrigerant: -101Kpa~4Kpa

Accuracy: $\pm 0.5\%$

Overload: 100Bar,10Mpa

Resolution: 1Kpa

Technical Specification :

Four Way Manifolds Digital



Laboratory Instrument India

