



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Standard Refrigeration System Training Device

**Product Code :**  
LBNY-0005-12100053



**Description :**

Standard Refrigeration System Training Device can be used to train students to master the theoretical knowledge and practice skills of the standard refrigeration system.

**Structure And Compositions:**

- Evaporator: tube air-cooled type, transparent chamber, and damper
- Compressor: 1/2HP, single phase 220V, control box
- Condenser: Air-cooled type, single phase 220V, include transparent glass tube
- Accumulator: 1/2HP, transparent glass type
- Filter dryer: 3/8 nut clamp type or welding type
- Sight glass: Welding type or nut clamp type
- Manometer: High pressure, low pressure
- Expansion valve: manual expansion valve
- Electronic valve: 3/8 nut clamp type or welding type.

**Technical Specification :**

Standard Refrigeration System Training Device

---

# Laboratory Instrument India

