



Email : sales@laboratoryinstrumentindia.com

Product Name :
Handheld LCR Meter Electronic Engineering Lab

Product Code :
LBNY-0002-620047



Description :

Handheld LCR Meter Electronic Engineering Lab

Technical Specification :

Handheld LCR Meter

26-segment analog bar graph; Series/parallel mode (selectable);

Open/short circuit calibration compensation & Smart selection mode for fast selection of components;

Mini-USB transmission to PC;

Loss factor (D), quality factor (Q), phase angle (?), equivalent series resistance (ESR), equivalent parallel resistance (EPR) and DCR DC resistance measurements.

Inductance (H): 20uH?2000H

Capacitance (F): 200pF?20mF

Resistance (?): 20??200M?

Handheld LCR Meter Electronic Engineering Lab, Handheld LCR Meter Electronic Engineering Lab Equipments, Handheld LCR Meter Electronic Engineering Lab Tools, Handheld LCR Meter Electronic Engineering Lab Tool Kits, Handheld LCR Meter Electronic Engineering Lab Manufacturers, Handheld LCR Meter Electronic Engineering Lab Suppliers from India, China, Kenya

Laboratory Instrument India

