



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Car Lucid Whole Model

**Product Code :**  
LBNY-0005-12100066



**Description :**

The interior mechanical structure and function of car is clearly understood under its demonstration and operation to make students and learners raise their initial knowledge of the Audi car's structure and its working principle. Car Lucid Whole Model is professional teaching equipment for colleges, secondary vocational and technical schools, car driving training and car repairing field.

**Technical Specification :**

Voltage: 220V single-phase; max power: 200W  
Overall size: Length 350cm \* Width 160cm \* Height \* 130cm  
Total net weight: 100kg.



**Laboratory Instrument India**