



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
DONGFENG Lucid Whole Model

**Product Code :**  
LBNY-0005-12100070



**Description :**

The interior mechanical structure and function of truck is clearly understood under its demonstration and operation to make students and learners raise their initial knowledge of the truck's structure and its working principle. DONGFENG Lucid Whole Model is professional teaching equipment for colleges, vocational schools, truck driving training and vehicle repairing field.

**Technical Specification :**

Voltage: 220V single-phase; max power: 200W  
Overall size: Length 320cm \* Width 125cm \* Height 130cm (60% real size)  
Total net weight: 100Kg.



**Laboratory Instrument India**