



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Transmission Model

**Product Code :**  
LBNY-0005-12100077



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The interior mechanical structure and function of car transmission is clearly understood through its demonstration and operation. Transmission Model is ideal teaching equipment for colleges, secondary vocational and technical schools, car driving training and vehicle repairing field.

**Technical Specification :**

Overall size: Length 50cm\*Width 30cm\*Height 25cm  
Total net weight: 3Kg



**Laboratory Instrument India**