



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Transmission Module Model

**Product Code :**  
LBNY-0005-12100084



**Description :**

The interior mechanical structure and function of gearbox is clearly understood under its demonstration and operation. Transparent Transmission educational training equipment is perfect teaching equipment for colleges, secondary vocational and technical schools, truck driving training and vehicle repairing field.

**Technical Specification :**

Overall size: Length 65cm \* Width 28cm \* Height 60cm  
Total net weight: 10kg.



**Laboratory Instrument India**