



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Truck Lucid Whole Model

**Product Code :**  
LBNY-0005-12100087



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The interior mechanical structure and function of truck is clearly understood under its demonstration and operation to make students and learners raise their initial knowledge of the truck's structure and its working principle. Transparent Truck educational equipment is professional teaching equipment for colleges, vocational schools, Truck driving training and vehicle repairing field.

**Technical Specification :**

Voltage: 220V single-phase; max power: 200W  
Overall size: Length 400cm \* Width 140cm \* Height \* 150cm  
Total net weight: 100kg.



**Laboratory Instrument India**