



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Magnetic Electric Starter Model

**Product Code :**  
LBNY-0005-121000106



**Description :**

The interior mechanical structure and function of magnetic starter is clearly understood under its demonstration and operation. Transparent Magnetic Electric Starter Model educational equipment is professional teaching equipment for colleges, secondary vocational and technical schools, car driving training and vehicle repairing field.

**Technical Specification :**

Overall size: Length 42cm \* Width 21cm \* Height \* 30cm  
Total net weight: 3kg



**Laboratory Instrument India**