



Email : sales@laboratoryinstrumentindia.com

Product Name :
Central Lock And Anti Theft System Teaching Board

Product Code :
LBNY-0005-121000111



Description :

It can fully display the composition and work process of the car alarm and central locking system. The teaching board contains car alarm and central locking system.

Features:

The panel is an aluminum panel. The color schematics is very convenient for students to analyses and to understand the work process.

It has the faults setting function: After setting the fault, the students can use apparatus to find the fault and do the fault troubleshooting exercises.

It is based on the real anti-theft device and central locking. The working principle and process of the anti-theft device and the central locking system can be demonstrated.

The electrical signals can be measured, such as resistance, voltage, current, frequency, etc of the electrical components.

Structure and composition:

Original car central locking ECU, glass lifting computer, glass lifter, four-door locks machine, tail cover lock, door lock switch, relay, data cable, high-capacity battery or 12V switch power, failure settings and troubleshooting system, color circuit schematics panel handling by computer, control panel, mobile stand.

Technical Specification :

Central Lock And Anti Theft System Teaching Board



Laboratory Instrument India