




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Integrated Ignition System Teaching Board

**Product Code :**  
LBNY-0005-121000113



#### **Description :**

The composition and the working process of various ignition system can be demonstrated. The device is based on the multiple ignition system.

#### **Structure and composition:**

The ignition system training device contains total six systems be below:

Magnetic induction ignition system,

Photoelectric induction ignition system

Platinum contact ignition system,

Holzer electronic ignition system,

Independent type ignition system without distributor

Group type ignition system without distributor

And include large capacity battery or 12V switch power, detection terminal, fault settings and troubleshooting system, mobile bench.

#### **Features:**

It has the specialized troubleshooting and set areas.

The color circuit diagrams are printed on the panel which is easy for the students to analyses and to understand the ignition system.

1-6 types ignition system can be demonstrated.

#### **Technical Specification :**

#### **Production technology:**

---

Power supply: DC12V, battery

Nylon casters with locking casters are assembled under the frame. it can be moved freely, and facilitate the teaching

Dimensions: about 1260x850x1860mm.



**Laboratory Instrument India**