




Email : sales@laboratoryinstrumentindia.com

Product Name :
Double Sided Basic Through Advanced Hydraulics Trainer

Product Code :
LBNY-0005-121000124



Description :

The Double Sided Basic Through Advanced Hydraulics Trainer has its unique ability to replicate a fully functional industrial or mobile hydraulic system, which is capable of operating with or without load.

Educational Advantages:

Allows students to:

Replicate a fully functional industrial or mobile hydraulic system, which is capable of operating with or without load

Operate the hydraulic system through complete load-cycles

Operate in fixed displacement, pressure compensated, and pressure/flow pump environments.

Technical Specification :

Includes:

Fixed displacement pump

Pressure-compensated pump

Hydraulic components

Diagnostic instruments

Load pump

Optional load-sensing and plug-n-play modules.

Laboratory Instrument India

