




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Double Sided Basic Through Advanced Hydraulics Trainer

**Product Code :**  
LBNY-0005-121000124



**Description :**

The Double Sided Basic Through Advanced Hydraulics Trainer has its unique ability to replicate a fully functional industrial or mobile hydraulic system, which is capable of operating with or without load.

**Educational Advantages:**

Allows students to:

Replicate a fully functional industrial or mobile hydraulic system, which is capable of operating with or without load

Operate the hydraulic system through complete load-cycles

Operate in fixed displacement, pressure compensated, and pressure/flow pump environments.

**Technical Specification :**

**Includes:**

Fixed displacement pump

Pressure-compensated pump

Hydraulic components

Diagnostic instruments

Load pump

Optional load-sensing and plug-n-play modules.

---

# Laboratory Instrument India

