



Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Sided Hydraulics Trainer With Touch Screen

Product Code :
LBNY-0005-121000125



Description :

This is an excellent feature for machine and equipment manufacturers who want to teach their personnel and clients how to setup and/or troubleshoot their own machines. The Double Sided Hydraulics Trainer With Touch Screen has the same impressive capability as the models for teaching a course in basic through advanced hydraulics. As an added option, you can program any machines hydraulic system into the PCs database. In addition to having all the practical hydraulics activities programmed into the onboard PC touch-screen, there is also an extraordinarily comprehensive diagnostics program.

Educational Advantages:

Allow students to:

Have all the practical hydraulics activities programmed into the onboard PC touch-screen

Teach a complete course in basic and advanced hydraulics

Program any machines hydraulic system into the PCs database

Troubleshoot hydraulic components and systems.

Technical Specification :

Include:

Optional load-sensing and plug-n-play modules

Panel PC with 15 color touch-screen.



Laboratory Instrument India