



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Front Auto Body Trainer Ford Focus

**Product Code :**  
LBNY-0005-121000134



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

This non-foldable trainer has been designed for auto body repair and refinishing classes. The Front Auto Body Trainer includes the hood and both front fenders.

**Main Features:**

Used for auto body repair and refinishing classes  
Includes hood and both front fenders.

**Technical Specification :**

**Educational Advantages:**

Takes up much less storage space than a real vehicle  
Allows real world repair and refinishing skill building without the requirement of having a complete vehicle in the shop.



**Laboratory Instrument India**