



Email : sales@laboratoryinstrumentindia.com

Product Name :
Front Auto Body Trainer Ford Focus

Product Code :
LBNY-0005-121000134



www.laboratoryinstrumentindia.com

Description :

This non-foldable trainer has been designed for auto body repair and refinishing classes. The Front Auto Body Trainer includes the hood and both front fenders.

Main Features:

Used for auto body repair and refinishing classes
Includes hood and both front fenders.

Technical Specification :

Educational Advantages:

Takes up much less storage space than a real vehicle
Allows real world repair and refinishing skill building without the requirement of having a complete vehicle in the shop.



Laboratory Instrument India