



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dual Door Trainer Non Powered

Product Code :
LBNY-0005-121000142



www.laboratoryinstrumentindia.com

Description :

Dual Door Trainer Non Powered trainer is designed for auto body and refinishing classes. The Dual Door Trainer includes the front doors and mirrors.

Main Features:

Used for auto body repair and refinishing classes
Includes doors and mirrors - non powered
Comes with interior door panels and hardware.

Technical Specification :

Educational Advantages:

Provides additional instructional experiences for door trim, door locks, window glass repair, door latch service, etc.
Provides real world learning experiences and skill development without the requirement of having a complete vehicle.
Takes up much less floor space than a real vehicle.



Laboratory Instrument India

