



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Dual Door Trainer Non Powered

**Product Code :**  
LBNY-0005-121000142



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Dual Door Trainer Non Powered trainer is designed for auto body and refinishing classes. The Dual Door Trainer includes the front doors and mirrors.

**Main Features:**

Used for auto body repair and refinishing classes  
Includes doors and mirrors - non powered  
Comes with interior door panels and hardware.

**Technical Specification :**

**Educational Advantages:**

Provides additional instructional experiences for door trim, door locks, window glass repair, door latch service, etc.  
Provides real world learning experiences and skill development without the requirement of having a complete vehicle.  
Takes up much less floor space than a real vehicle.



**Laboratory Instrument India**

