



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cutaway Over Center Variable Displacement Pump

**Product Code :**  
LBNY-0005-121000143



www.laboratoryinstrumentindia.com

**Description :**

The Cutaway Over Center Variable Displacement Pump trainer is an application that can be chosen to demonstrate the operation of a hydraulic over-center pump.

**Educational Advantages:**

The pump housing is cutaway to expose the internal components  
The barrel is cutaway to expose the pistons  
The rear cover is cutaway to expose the pressure relief valve  
The pump is lightweight and easy to handle  
The swash plate moves  
Large knob on output shaft makes it easy to rotate  
The charge pump housing is cutaway to expose the charge pump and the charge pressure relief valve  
Includes stand.

**Technical Specification :**

Cutaway Over Center Variable Displacement Pump



**Laboratory Instrument India**

