



Email : sales@laboratoryinstrumentindia.com

Product Name :
Squirrel Cage Induction Machine

Product Code :
LBNY-0005-121000146



www.laboratoryinstrumentindia.com

Description :

This Squirrel Cage Induction Machine 60 HZ is a three-phase squirrel cage induction motor. The three-stator windings may be connected in wye or delta and be brought out to separate 4mm terminals on the housing. Specifications and/or diagrams are silk-screened on the back of a clear-plastic front panel for ease of identification. The B design of the machine makes it convenient for most applications. It can be used as either a motor or generator. The trainer is protected by a either a diode or thermal-magnetic circuit breaker. Included in the Rotating Machines Trainer for 1/3 hp - 250 W. Requires a 120V 60Hz circuit.

Educational Advantages:

The three-stator windings may be connected in wye or delta
Can be used as either a motor or generator
Included in the basic Rotating Machines Trainer for 1/3 hp - 250 W.

Technical Specification :

Squirrel Cage Induction Machine



Laboratory Instrument India

