



Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Phase AC DC Supply 60 HZ

Product Code :
LBNY-0005-121000147



www.laboratoryinstrumentindia.com

Description :

The module is housed into a half size cabinet which is totally enclosed, portable and modular. Connections to the unit are made through 4 mm standard color-coded, five-way binding posts. The trainer is a Single Phase Power Supply unit used for various experiments requiring variable or fixed AC and DC voltages. Specifications and/or diagrams are silk-screened on the back of a clear-plastic front panel for ease of identification. This permits the utilization of banana-type jumper leads between the single phase power supply and other components in the circuit. The trainer is protected by a either a diode or thermal-magnetic circuit breaker. Requires a 120V 60Hz circuit.

Educational Advantages:

Fixed and variable single phase AC output
Single phase input
Variable DC output.

Technical Specification :

Single Phase AC DC Supply 60 HZ



Laboratory Instrument India

