



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fender Support

Product Code :
LBNY-0005-121000149



www.laboratoryinstrumentindia.com

Description :

The Fender Support is used for holding fenders during working, buffing and refinishing tasks in auto body classes.

Main Features:

Comes as a set of two
Holds any type of fender solidly in place on workbenches and tables
Used for fender working, buffing and refinishing

Educational Advantages:

Portable and stores easily when not in use
Solves the common problem of providing a secure and steady method of holding fenders for body repair procedures.

Technical Specification :

Fender Support



Laboratory Instrument India

