



Email : sales@laboratoryinstrumentindia.com

Product Name :
Communication Trainer Electronic Engineering Lab

Product Code :
LBNY-0002-630004



Description :

Communication Trainer Electronic Engineering Lab

Technical Specification :

Communication Trainer
Function Generator
Complete dual channels function/arbitrary waveform generator
Channel independence, coupling, track working mode
200 MSa/s sampling rate and 14-bit vertical resolution per channel
Output of 6 standard waveforms, built-in 50 kinds of arbitrary waveform
Multi modulation function: AM, DSSC - AM, FM, PM, PWM, FSK, ASK,
BPSK and logarithm/linear sweep 0-250 V, 2A, with Voltmeter
DC power supply panel, adjustable 0-30 V, can supply power 2A
Dual Universal out let power supply panel 2P+PE size 220 Volt Use with a current of 10 amps, 1 set
Function Generator panel size 0.5Hz-2.5MHz can generate sine, triangle.

Communication Trainer Electronic Engineering Lab, Communication Trainer Electronic Engineering Lab Equipments, Communication Trainer Electronic Engineering Lab Tools, Communication Trainer Electronic Engineering Lab Tool Kits, Communication Trainer Electronic Engineering Lab Manufacturers, Communication Trainer Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India