



Email : sales@laboratoryinstrumentindia.com

Product Name :
Side Auto Body Trainer

Product Code :
LBNY-0005-121000151



www.laboratoryinstrumentindia.com

Description :

This trainer has been designed for auto body and refinishing classes and can be made on any model of your choice. The Side Auto Body Trainer includes a front fender, door(s) and rear quarter panel.

Main Features of Side Auto Body Trainer
Used for auto body repair and refinishing classes
Vehicle of your choice
Complete side of a vehicle body
Includes front fender, door(s) and rear quarter panel.

Technical Specification :

Educational Advantages:

Provides additional instructional experiences for door trim, door locks, window regulator/glass repair, door latch service, etc.
Provides real world learning experiences and skill development without the requirement of having a complete vehicle.
Takes up much less floor space than a real vehicle.

Laboratory Instrument India

