



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Side Auto Body Trainer

**Product Code :**  
LBNY-0005-121000151



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

This trainer has been designed for auto body and refinishing classes and can be made on any model of your choice. The Side Auto Body Trainer includes a front fender, door(s) and rear quarter panel.

Main Features of Side Auto Body Trainer  
Used for auto body repair and refinishing classes  
Vehicle of your choice  
Complete side of a vehicle body  
Includes front fender, door(s) and rear quarter panel.

**Technical Specification :**

**Educational Advantages:**

Provides additional instructional experiences for door trim, door locks, window regulator/glass repair, door latch service, etc.  
Provides real world learning experiences and skill development without the requirement of having a complete vehicle.  
Takes up much less floor space than a real vehicle.

---

# Laboratory Instrument India

