



Email : sales@laboratoryinstrumentindia.com

Product Name :
Contactless Torque Meter Refurbishing

Product Code :
LBNY-0005-121000154



www.laboratoryinstrumentindia.com

Description :

Refurbishing your old-version Contactless Torque Meter is a cost-effective alternative to replacing costly equipment.

Refurbishing Work

Replacement Of Torque Sensor (Transducer)
Adding of an analog output mounted on a BNC plug
Adding of two DEL indicators for clockwise/counterclockwise
Installation of electronic contactless transducer (50/0-30 Nm)
Adding of an electronic board for digital torque readout

Complete Console Disassembling
Housing mechanical repair if required
Sand blasting and painting

Complete Testing And Verifying
Inspection report and quality control
Electrical and mechanical testing.

Technical Specification :

Contactless Torque Meter Refurbishing



Laboratory Instrument India