



Email : sales@laboratoryinstrumentindia.com

Product Name :
Contactless Torque Meter Refurbishing

Product Code :
LBNY-0005-121000154



www.laboratoryinstrumentindia.com

Description :

Refurbishing your old-version Contactless Torque Meter is a cost-effective alternative to replacing costly equipment.

Refurbishing Work

- Replacement Of Torque Sensor (Transducer)
- Adding of an analog output mounted on a BNC plug
- Adding of two DEL indicators for clockwise/counterclockwise
- Installation of electronic contactless transducer (50/0-30 Nm)
- Adding of an electronic board for digital torque readout

- Complete Console Disassembling
- Housing mechanical repair if required
- Sand blasting and painting

- Complete Testing And Verifying
- Inspection report and quality control
- Electrical and mechanical testing.

Technical Specification :

Contactless Torque Meter Refurbishing



Laboratory Instrument India