



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrical Motor Circuit Troubleshooting 60 HZ

Product Code :
LBNY-0005-121000158



Description :

The training panel consists of: All industrial components required to complete lab experiments; One power cord with plug for 120/208 V AC - 3 phases; One manual fault box with forty-four various faults; One user manual with panel wiring diagrams; Two 1/3 HP (250 W) electric motors including one squirrel cage motor and one capacitor start motor. The Electrical Motor Circuit Troubleshooting 60 HZ is designed for the electric motor circuit troubleshooting tasks.

Main Features:

18 Gauge steel panel
Lockable manual fault box fault simulator including 44 operating faults.
Operating manual with wiring diagrams.
Two base-fixed stands for safety installation on working tables.

Technical Specification :

Electrical Motor Circuit Troubleshooting 60 HZ



Laboratory Instrument India

