



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
AC Ammeters 50/60 HZ

**Product Code :**  
LBNY-0005-121000159



www.laboratoryinstrumentindia.com

**Description :**

The AC Ammeters 50/60 HZ trainer includes three analog ammeters for reading alternative current. The trainer is housed into a totally enclosed, stackable and portable, full-size modular cabinet. Protected by diodes. Schematic diagrams and specifications are silk-screened onto the back of the front panel of each unit. Interconnections between components and connections are made through standard color coded safety sockets with banana plug jumper leads. Each ammeter has four reading scales that are adapted to measure the current in experiments with rotary machines, electrical circuits or motor controllers. Included in the Rotating Machines Trainers for 1/3 hp - 250 W.

**Educational Advantages:**

Each ammeter has four reading scales that are adapted to measure the current in experiments with rotary machines, electrical circuits or motor controllers  
Includes three analog ammeters for reading alternative current.

**Technical Specification :**

Taut-band type: left zero  
Diode or thermal protection  
AC ammeters (3): 0-1 / 3 / 10 / 30 A AC  
Accuracy:  $\pm 2\%$ .

---

# Laboratory Instrument India

