



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Rockwell Hardness Tester Electronic Engineering Lab

**Product Code :**  
LBNY-0002-640009



**Description :**

Rockwell Hardness Tester Electronic Engineering Lab

**Technical Specification :**

Hardness scale HRA,HRB,HRC,HRD,HRF,HRG  
Preliminary test force: 98N (10KG)  
Test force: 588n(60kg) 980n (100kg) 1471n (150kg)  
Display: Analogue  
Stage elevation: Manual  
Load control: Manual  
Reading: 0.5HR  
Work place height: 170mm  
Workpiece depth: 135mm (from center of indenter )  
Dimension: 466x238x630mm  
Weight: 70kg

Rockwell Hardness Tester Electronic Engineering Lab, Rockwell Hardness Tester Electronic Engineering Lab Equipments, Rockwell Hardness Tester Electronic Engineering Lab Tools, Rockwell Hardness Tester Electronic Engineering Lab Tool Kits, Rockwell Hardness Tester Electronic Engineering Lab Manufacturers, Rockwell Hardness Tester Electronic Engineering Lab Suppliers from India, China, Kenya



**Laboratory Instrument India**