



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rockwell Hardness Tester Electronic Engineering Lab

Product Code :
LBNY-0002-640009



Description :

Rockwell Hardness Tester Electronic Engineering Lab

Technical Specification :

Hardness scale HRA,HRB,HRC,HRD,HRF,HRG
Preliminary test force: 98N (10KG)
Test force: 588n(60kg) 980n (100kg) 1471n (150kg)
Display: Analogue
Stage elevation: Manual
Load control: Manual
Reading: 0.5HR
Work place height: 170mm
Workpiece depth: 135mm (from center of indenter)
Dimension: 466x238x630mm
Weight: 70kg

Rockwell Hardness Tester Electronic Engineering Lab, Rockwell Hardness Tester Electronic Engineering Lab Equipments, Rockwell Hardness Tester Electronic Engineering Lab Tools, Rockwell Hardness Tester Electronic Engineering Lab Tool Kits, Rockwell Hardness Tester Electronic Engineering Lab Manufacturers, Rockwell Hardness Tester Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India