



Email : sales@laboratoryinstrumentindia.com

Product Name :
Squirrel Cage Machine Refurbishing

Product Code :
LBNY-0005-121000165



www.laboratoryinstrumentindia.com

Description :

Refurbishing your old-version Unitec Squirrel Cage Induction Machine is a cost-effective alternative to replacing costly equipment.

Refurbishing Work
Console Mechanical Reassembling
Console reassembling
Working table replacement
Clear rear panel changing
Repair and/or replacement of casters and couplings
Assembling and mechanical adjustment of motor

Complete Console Disassembling
Mechanical repair of housing if required
Sand blasting and painting

Complete Testing And Verifying
Coupling verifying on both sides of machine
Inspection report and quality control
Electrical testing w/o load

Electrical Reassembling And Rewiring
Front panel replacement
Rewiring

Electrical part replacement

Complete Motor Disassembling

Motor inspection

Repair and replacement of brush assembly if required

External painting

Bearing replacement.

Technical Specification :

Squirrel Cage Machine Refurbishing



Laboratory Instrument India