



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wound Rotor Induction Machine Refurbishing

Product Code :
LBNY-0005-121000171



www.laboratoryinstrumentindia.com

Description :

Refurbishing your old-version Unitec Wound Rotor Induction Machine is a cost-effective alternative to replacing costly equipment.

Refurbishing Work

Complete Testing And Verifying
Coupling verifying on both sides of machine
Inspection report and quality control
Electrical testing w/o load

Complete Console Disassembling
Mechanical repair of housing if required
Sand blasting and painting

Complete Motor Disassembling
Motor inspection
Repair and replacement of brush assembly if required
External painting
Bearing replacement

Electrical Reassembling And Rewiring
Front panel replacement
Electrical part replacement
Rewiring.

Console Mechanical Reassembling
Console reassembling
Clear rear panel changing
Assembling and mechanical adjustment of motor
Repair and/or replacement of casters and couplings
Working table replacement.

Technical Specification :

Wound Rotor Induction Machine Refurbishing



Laboratory Instrument India