



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mechanical Load Mounting Base

**Product Code :**  
LBNY-0005-121000174



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It is used to hold the units included in the Mechanical Load Assembly. Can also be inserted in the machine coupling bases. The Mechanical Load Mounting Base is made from rigid black laminate and is equipped with carefully spaced mounting holes to ensure safe connections.

**Educational Advantages:**

Equipped with wing screw anchors to safely bolt the units in the experiments  
Rigid laminated wooden base.

**Technical Specification :**

Material: Laminated wood, black  
Number of screw anchors: 4  
Screw length: 5 cm (2")  
Dimensions: 25 cm x 31.25 cm x 20 cm (10 " x 12 ½" x 8")  
Weight: 1.36 kg (3 lbs).



**Laboratory Instrument India**

