



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
AC Voltmeters 60 HZ

**Product Code :**  
LBNY-0005-121000176



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The unit is housed in a totally enclosed, stackable and portable, half-size modular cabinet. Protected by thermal magnetic breakers. Schematic diagrams and specifications are silk-screened onto the back of the front panel of each unit. The AC Voltmeters 60 HZ trainer includes three analog voltmeters for reading AC voltages. Interconnections between components and connections are made through standard color coded safety sockets with banana plug jumper leads. Each voltmeter has two reading scales adapted to measure the voltage in experiments with rotary machines, electrical circuits, or motor controllers. Included in the Rotating Machines Trainer for 1/3 hp - 250 W.

**Educational Advantages:**

Each voltmeter has two reading scales adapted to measure the voltage in experiments with rotative machines, electrical circuits, or motor controllers  
Includes three analog voltmeters for reading AC voltages.

**Technical Specification :**

Diode or thermal protection  
AC voltmeters (3): 0-150 / 300 V AC  
Accuracy:  $\pm 2\%$ .

---

# Laboratory Instrument India

