



Email : sales@laboratoryinstrumentindia.com

Product Name :
Synchronous Machine Refurbishing

Product Code :
LBNY-0005-121000184



www.laboratoryinstrumentindia.com

Description :

Refurbishing Work
Console Mechanical Reassembling
Console reassembling
Clear rear panel changing
Assembling and mechanical adjustment of motor
Repair and/or replacement of casters and couplings
Working table replacement

Complete Console Disassembling
Mechanical repair of housing if required
Sand blasting and painting

Complete Testing And Verifying
Coupling verifying on both sides of machine
Electrical testing w/o load
Inspection report and quality control.

Complete Motor Disassembling
Motor inspection
Repair and replacement of brush assembly if required
External painting
Bearing replacement

Electrical Reassembling And Rewiring
Electrical part replacement
Front panel replacement
Rewiring.

Technical Specification :

Synchronous Machine Refurbishing



Laboratory Instrument India